**INSTRUCTIONS FOR USE**

CLARAVUE LONG TERM  
PATIENT MONITORING

**CE** EC ACCORDING TO  
EU MDR 2017/745

**i** READ INSTRUCTIONS  
FOR USE

**M** MANUFACTURER

**REF** REFERENCE

**LOT** LOT CODE

**Use By** EXPIRY

**TEMPERATURE LIMITS**

**KEEP AWAY FROM  
DIRECT SUNLIGHT**

**PROTECT FROM  
MOISTURE**

**DO NOT USE IF  
PACKAGING IS DAMAGED**

**MR** MR UNSAFE

**DURATION OF USE**

**SINGLE PATIENT USE**

**MD** MEDICAL DEVICE

**UDI** UNIQUE DEVICE  
IDENTIFICATION

**QTY** QUANTITY

**PHT** PHT FREE

**NOT MADE WITH  
NATURAL RUBBER  
LATEX**

**NOT MADE WITH BPA**

**ROHS** THE PRODUCT IS  
ROHS COMPLIANT

**DO NOT THROW  
THE PRODUCT INTO  
COMMON WASTE.  
RESPECT THE  
NATIONAL  
REGULATION.**

**Rx ONLY** CAUTION: U.S.  
FEDERAL LAW  
RESTRICTS THIS  
DEVICE TO SALE BY  
OR ON THE ORDER  
OF A PHYSICIAN

**X-RAY TRANSLUCENT**

HS.NisshaMedical.com/IFU

**CLARAVUE LONG TERM PATIENT MONITORING**

**I. DESCRIPTION AND SPECIFICATION**

**INTENDED USE**

**CLARAVUE® LPM** are prewired electrodes intended for use in electrocardiographic procedures as passive sensors for the derivation of electrical signals from the heart in a hospital setting by qualified medical personnel. The prewired electrodes are pre-gelled, single patient use, non-sterile, and are to be applied on clean intact skin. Intended patient groups are adults.

**IDENTIFICATION**

Part number	Description	Population	X-Ray Compatibility
50700	CLARAVUE LPM 3 RT (SIDE)	Adults	Yes
50701	CLARAVUE LPM 5 RT (SIDE)	Adults	Yes
50720	CLARAVUE LPM 3 RT (TOP)	Adults	Yes
50721	CLARAVUE LPM 5 RT (TOP)	Adults	Yes

**CLARAVUE® LPM** product family meets the requirements of the standard, ECG TRUNK CABLES and PATIENT LEADWIRES (ANSI/AAMI EC53). **CLARAVUE® LPM** is intended for SINGLE-PATIENT USE.

**PERFORMANCE CHARACTERISTICS**

<b>Gel type</b>	Adhesive hydrogel
<b>Sensor</b>	Ag/AgCl
<b>Adapter/Connector</b>	Integrated leadwire
<b>Compatibility</b>	IP Cables & Leadwires product range
<b>Wear time</b>	Up to 120 hours
<b>Shelf-life</b>	24 months

**II. GENERAL WARNING AND SAFETY PRECAUTIONS  
CONTRAINDICATIONS FOR USE**

There are no specific contraindications known for using the product.

**Warnings and safety precautions**

- **CLARAVUE® LPM** is for single-patient-use only and must not be used on multiple patients. Use on other patients can cause patient infection or cross-infection.
- **CLARAVUE® LPM** must be used on intact (non-injured) skin only.
- MR Unsafe - These **CLARAVUE® LPM** should not be operated in magnetic resonance imaging environments
- These **CLARAVUE® LPM** should not be operated with monophasic defibrillators
- This product can only be used by trained personnel familiar with the proper placement and use of the electrodes.
- Make sure to carefully read this instruction for use, as well as the ones of the other medical devices that may be used along.
- Make sure to check the integrity of the packaging before use. Look for tears on the pouch or the sealing. Discard the product if there are any marks of damage.
- Do not use a product found in packaging which expiry date has passed.
- Do not modify the product.
- Do not soak, rinse, or sterilize this device as it may cause malfunction of the device
- U.S. federal law restricts this device to sale by or on the order of a physician

Report any serious incident that occurred in relation to the device to the manufacturer at **FRquality@nisshamedical.com** and to the competent authority of the Member State where the hospital is established.

**II. INSTALLATION / USE / CLEANING / MAINTENANCE / DISPOSAL****SELF HANDLING AND PREPARATION:**

- Remove excess hair from the electrode sites to ensure maximal contact
- Clean the electrode sites with mild soap and water or a non-alcoholic wipe
- Make sure that the skin is completely dry before the application of the electrodes.
- If necessary, perform a light exfoliation of the skin gently to improve the adherence of the product.

**APPLICATION/ USE:**

- Check the product for physical damage before use.
- The product should be replaced if the electrodes no longer stick firmly to the skin.

1. Open the single **CLARAVUE® LPM** patient kit only immediately before use. Remove inner bag. Unwind sets from the holder card.
2. Check if the electrodes are undamaged. In case of any damage, discard the product.
3. Gently remove electrode sets from the release liner by peeling away from the thumb tab.
4. Place the electrodes on the patient.
5. Connect the set of Leadwires to the appropriate trunk cable from **IP Cables & Leadwires** product range.
6. Start ECG recording. If the ECG trace is not acceptable, check that the electrodes are positioned correctly and well attached to the skin. If the trace is still not acceptable, discard the product and place a new one.

**REMOVAL:**

- Uncareful removal of electrodes can cause skin irritation.

1. Unplug the set of leadwires from the trunk cable.
2. Carefully remove the electrodes from the patient's skin. Keep the electrode close to the skin surface as you pull it back and support the skin with a finger immediately adjacent to the adhesive being removed.
3. If the electrodes are difficult to remove, use water to dissolve the adhesive.
4. Gel residue can be removed using paper towels or water.
5. To minimize the risk of skin irritation, find a new area of skin on which to place the electrodes when repositioning the product.
6. Throw away the product after use (see "Disposal").

**PACKAGING/STORAGE:**

- The **CLARAVUE® LPM** should be stored in the original packaging, in temperature ranging from 15 to 35° C, protected from the sunlight, dampening, and dirt. When opened, the product must be used immediately.

**MANUAL CLEANING / DISINFECTING:**

- The **CLARAVUE® LPM** must not be reprocessed. Please discard after use.

**MAINTENANCE:**

- The **CLARAVUE® LPM** must not be repaired. Please discard any defective product.

**DISPOSAL:**

- Discard the product after duration of use on one patient, or immediately discard if any damage or deterioration is observed.
- The **CLARAVUE® LPM** should be disposed of according to hospital waste policies and procedures.

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