



TenderTrode™

Tender Care for Tender Patients

FEATURES



INFANT BREATHABLE ELECTRODES

Designed for patient comfort. Gentle adhesive and low chloride hydrogel provide a conforming and non-irritating electrode. Excellent for use in APNEA monitoring and all other pediatric applications.



SOFT CLOTH ELECTRODES

Adhesive solid gel electrodes with soft cloth material are available in various sizes to suit your needs.



FOAM PEDIATRIC ELECTRODES

Our foam pediatric electrodes are available in 1" diameter to allow for easy placement and features solid gel that is gentle on the skin.

For use in NICU, Diagnostics, Apnea and PICU.

PACKAGING



RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes



CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | HS.NisshaMedical.com

Patient Monitoring & Medical Surgical Consumables

Latex-free



Solid Gel



Wet Gel



TenderTrode™

Tender Care for Tender Patients

Vermed TenderTrode™ electrodes from NMT - HS represent a complete selection of electrodes for infant and pediatric use. Each electrode is carefully designed with gentle, hypoallergenic adhesives and soft conforming materials to ensure comfort and effectiveness on young patients.

A10001 - 1.625" / 41.275mm - ROUND TAPE ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10001	3	3	n/a	600		
A10000-60	3	60	n/a	600		72 Hours

A10040 - 0.875" / 22.225mm - SQUARE CLOTH ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10040-60	3	60	n/a	600		

A10041 - 0.875" x 1.3125" / 22.225mm x 33.34mm - RECT. CLOTH ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10041-60	3	60	n/a	600		

A10042 - 1.375" / 34.925mm - SQUARE CLOTH ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10042-60	3	60	n/a	600		

P7-3 - 0.875" / 22.225mm - TEARDROP FOAM ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A1P7-3-60	3	3	108	1080		

A10027 - 0.5" x 0.9 / 12.7mm x 22.86mm - AG/AGCL TinyTab ELECTRODE



Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10027	10	100	n/a	1000		