



## TenderTrode™

Pre-Wired Electrodes & Thermal Probe Covers

### FEATURES

#### PRE-WIRED ELECTRODES



We offer a complete line of pre-wired electrodes for use in PICU/NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.


#### THERMAL PROBE COVERS



Our thermal probe covers use a delicate adhesive designed for tender infant skin. The highly reflective surface deflects heat from radiant warmers.

A complete line of pre-wired electrodes for use in the PICU/NICU and other neonatal applications.

### PACKAGING



#### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

#### CASE

FDA / GS1 Compliant Labels



# TenderTrode™

Pre-Wired Electrodes & Thermal Probe Covers

Nissha Medical Technologies offers a complete line of pre-wired electrodes for use in the PICU/ NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.

## A10057-S - 0.875" / 22.225mm - SQUARE PRE-WIRED FOAM ELECTRODE



Part #	Product #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6502910V	A10057-S	3	3	300		 72 Hours

## A10061-SRT - 0.875" / 22.225mm - SQUARE PRE-WIRED CLOTH ELECTRODE - RT



Part #	Product #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6502911V	A10061-SRT	3	3	300		 72 Hours

## A10010-S - 0.875" X 1.875" / 22.225mm x 47.625mm - RECT. PRE-WIRED CLOTH ELECTRODE



Part #	Product #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6502912V	A10010-S	3	3	300		 72 Hours

## A10060-SRT - 0.875" X 1.875" / 22.225mm x 47.625mm - RECT. PRE-WIRED CLOTH ELECTRODE - RT



Part #	Product #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6502914V	A10060-SRT	3	3	300		 72 Hours

## A10008 - THERMAL PROBE COVER



Part #	Product #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
6500507V	A10008	1	50	500		 72 Hours