



Defibrillator Compatibility

Match Defibrillator Brands to Vermed Part Numbers

Wet Gel MR Conditional
 

| BRAND OF DEFIBRILLATOR | MODEL OF DEFIBRILLATOR | VERMED COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS | |
|------------------------|---|--|--|
| Philips Medical | <ul style="list-style-type: none"> • XLT • HeartStart XL • HP Codemaster XL (Plug Style Connector) | <ul style="list-style-type: none"> • MRX • MRX ALS | <ul style="list-style-type: none"> • D311 • D309 |
| | <ul style="list-style-type: none"> • HeartStart FR2+ | <ul style="list-style-type: none"> • HeartStart Forerunner | <ul style="list-style-type: none"> • D309 • D311 |
| | <ul style="list-style-type: none"> • HeartStart FR3 | <ul style="list-style-type: none"> • HeartStart FRx | NOT COMPATIBLE |
| Physio-Control | <ul style="list-style-type: none"> • LP-9P • LP-10C • LP-12 | <ul style="list-style-type: none"> • LP-15 • LP-20 | <ul style="list-style-type: none"> • D304 • D313 |
| | <ul style="list-style-type: none"> • LP-500 • LP-1000 | <ul style="list-style-type: none"> • LP-500 DPS | <ul style="list-style-type: none"> • D304 • D313 |
| Zoll Medical | <ul style="list-style-type: none"> • M-Series • R-Series | <ul style="list-style-type: none"> • E-Series • X-Series | <ul style="list-style-type: none"> • D303 • D312 |
| | <ul style="list-style-type: none"> • Zoll AED PLUS | | NOT COMPATIBLE |
| Cardiac Science | <ul style="list-style-type: none"> • PowerHeart AED G3 Plus • PowerHeart AED G3 Pro | | <ul style="list-style-type: none"> • D326 |

Please contact Vermed directly for any questions regarding compatibility.