


**STANDARD**  
ELECTRODES & LEADWIRES

CLARMUE®

**ECG ARTIFACT**  
Heavy leadwires and snaps cause increased motion artifact.

**REDUCED ECG ARTIFACT**  
Lightweight leadwires and pre-wired electrodes result in clearer traces.




**STANDARD ELECTRODES**  
Ten individual electrodes

**TIME TO APPLY**  
90 seconds


**HEAVY LEADWIRES**  
0.405 kg

**LEADWIRE CONNECTION CABLE**  
Ten individual leadwires require individual connection to trunk cable



**PRE-POSITIONED PAIRS**  
Eliminates the need to place each electrode individually

**TIME TO APPLY**  
45 seconds



**LEADWIRE CONNECTION CABLE**  
Patented electrode connector reduces connections to trunk cable

**LIGHTWEIGHT LEADWIRES**  
0.265 kg

STANDARD	COMPARISON	CLARMUE®
	Quality & Clear Trace	✓
0.405 kg	Lightweight System	0.265 kg
	Hygienic Disposable Leadwires	✓
	Reusable Trunk Cable	✓
	Radiolucent	✓
90 seconds	Quick and Easy to Apply	45 seconds
	Prewired Solution	✓





## CLARAVUE®

### Patient Monitoring Kit

Preferred by leading Cardiologists, Electrophysiologists, Cardiology Department Directors, and Cardiology Nurse Leaders, CLARAVUE® delivers the best-in-class, hygienic solution for extremely precise, radiolucent patient monitoring.

#### CLARAVUE® RADIOLUCENT DISPOSABLE PRE-WIRED ELECTRODES

Part #	Product #	Description	Per Kit	Kits/Box	Kits/Case	Solid Gel	Duration of Use
50401-US	32028776	Pediatric Patient Kit - 3 electrodes	3	40	480		 24 Hours
50600-US	32028777	Adult Patient Kit - 4 electrodes	4	40	480		
50601-US	32028778	Adult Patient Kit - 5 electrodes	5	30	360		
50602-US	32028779	Adult Patient Kit - 10 electrodes	10	25	300		
50603-US	32028780	Adult Patient Kit - 6 electrodes - V-Series	6	30	360		
50604-US	32028781	Adult Patient Kit - 5 electrodes - V-Series	5	30	360		

