

CLARAVUE®



HEARTSYNC
Defibrillation • Cardioversion • Pacing • ECG Monitoring




A COMPLETE MONITORING & DEFIBRILLATION SYSTEM FOR EP & CATH LABS

For the leading Cardiologists, Electrophysiologists, Cardiology Department Directors, and Cardiology Nurse Leaders, CLARAVUE® is the best-in-class, hygienic solution for extremely precise, radiolucent patient monitoring. *Please Note: CLARAVUE® is NOT defibrillation proof.*

Nissha Medical's HeartSync Defibrillation Electrodes are designed with your needs in mind. Featuring the longest leadwires, 36 month shelf life, and ease of placement, they will meet the needs of your clinical staff.

- HYGIENIC DISPOSABLE LEADWIRES**
Disposable Patient Kit ensures the most hygienic solution for infection control
- QUICK & EASY TO APPLY**
Patented quick connection point, repositionable & comfortable
- PREWIRED SOLUTION**
Eliminates the need to snap leadwires onto the individual electrodes
- FULLY RADIOLUCENT**
No image obstruction under x-ray or fluoroscopy
- QUALITY & CLEAR TRACE**
Prewired electrodes result in clearer trace
- LIGHTWEIGHT LEADWIRES**
0.265 Kg compared to 0.405 kg on standard leadwires

- 60 INCH LEAD WIRES**
LONGEST ON THE MARKET
Increased flexibility with equipment placement & clinician mobility
- 36 MONTH SHELF LIFE**
LONGEST ON THE MARKET
Extended shelf life reduces waste and saves money
- SMALLER PAD SIZE**
Allows for easier placement without crowding
- ADVANCED HYDROGEL**
Unique gel provides maximum adhesion
- LARGER CONDUCTOR PLATE**
Larger conductive area disperses energy more effectively
- COLOR CODED PACKAGING**
The correct product, for your device
- COLOR CODED PADS**
Anterior electrode is blue foam and the posterior electrode is white foam. This helps train staff to look for colors to aid in pad placement.
- COMPREHENSIVE OFFERING**
Radiolucent and Radiotransparent leads in/out and pediatric options available
- INSTRUCTIONS FOR USE**
User checklist, on the release liner. Liner is opaque and visible when dropped on floor
- LEADS OUT OPTION**
Allows pads to be pre-connected prior to opening pouch.



PHYSIO CONTROL

ZOLL MEDICAL

CARDIAC SCIENCE

ORDER INFORMATION

CLARAVUE®


RADIOLUCENT DISPOSABLE PREWIRED ELECTRODES



Product #	Description	Per Kit	Kits/Box	Kits/Case	Gel Type	Duration of Use
50401-US	Pediatric Patient Kit - 3 electrodes	3	40	480	Solid Gel	24 hours
50600-US	Adult Patient Kit - 4 electrodes	4	40	480	Solid Gel	24 hours
50601-US	Adult Patient Kit - 5 electrodes	5	30	360	Solid Gel	24 hours
50602-US	Adult Patient Kit - 10 electrodes	10	25	300	Solid Gel	24 hours
50603-US	Adult Patient Kit - 6 electrodes - V-Series	6	30	360	Solid Gel	24 hours
50604-US	Adult Patient Kit - 5 electrodes - V-Series	5	20	360	Solid Gel	24 hours



DEFIBRILLATION ELECTRODES



Product #	For Defibrillator Models	Leads	Pads & Leadwires
C100AC-PHYSIO	Lifepak 15, 20/20e, 1000	In	
C100LOAC-PHYSIO	Lifepak 15, 20/20e, 1000	Out	
T100AC-PHYSIO	Lifepak 15, 20/20e, 1000	In	
T100LOAC-PHYSIO	Lifepak 15, 20/20e, 1000	Out	
C100AC-ZOLL	R-Series, X-Series	In	
C100LOAC-ZOLL	R-Series, X-Series	Out	
T100AC-ZOLL	R-Series, X-Series	In	
T100LOAC-ZOLL	R-Series, X-Series	Out	
T100-CS	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	Can be used in place of Cardiac Science Part #9131-001

EP/CATHLAB



DEFIBRILLATION & CLARAVUE PREWIRED ELECTRODES

NISSHA
MEDICAL TECHNOLOGIES

03-13-01-066 (05/2023)