



# NeuroPlus<sup>™</sup>

Sleep Study Electrodes

## FEATURES:



### SOLID GEL

Specially formulated solid gel ensures full sensor contact and excellent readings.



### GENTLE ADHESIVE

Designed with a gentle adhesive gel for use in sensitive areas.



### REPOSITIONABLE

Specially formulated adhesive gel is repositionable to ensure proper placement.

For use in sleep studies, EMG, NCS and EP

## PACKAGING



### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

### CASE

FDA / GS1 Compliant Labels






# NeuroPlus™




Repositionable, Disposable Electrodes

NeuroPlus™ electrodes are smaller in size, allowing for easy placement on the face or other difficult locations. NeuroPlus™ features a cloth base or foam option for added comfort to the patient.




## A10040 - 0.875" / 22.225mm - SQUARE CLOTH ELECTRODE - REPOSITIONABLE

	Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
	6500060V	A10040-60	3	60	600		 48 Hours

## A10041 - 0.875" x 1.3125" / 22.225mm x 33.34mm - RECT. CLOTH ELECTRODE - REPOSITIONABLE

	Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
	6500061V	A10041-60	3	60	600		 48 Hours

## A10042 - 1.375" / 34.925mm - SQUARE CLOTH ELECTRODE - REPOSITIONABLE

	Product #	Ref #	Per Strip	Per Pouch	Per Case	Solid Gel	Duration of Use
	6500062V	A10042-60	3	60	600		 48 Hours