




VersaTrode™  
Adult General Purpose Electrodes

**FEATURES**



**WET & SOLID GEL ELECTRODES**

Wet gel provides quick reading and stabilization. Solid gel improves adhesion and requires minimal clean up.



**A10083 VERSATRODE™ REPOSITIONABLE**

MR Conditional certified carbon snap and unique solid gel allows for the electrode to be repositioned up to 5 times.

VersaTrode™ electrodes are the ideal choice for adult cardiac monitoring and diagnostic applications.

**PACKAGING**

**RESEALABLE BULK POUCH**

Ziplock pouch contains 30 or more electrodes



**CASE**

FDA / GS1 Compliant Labels



**HYGIENIC, DISPOSABLE SINGLE-USE PACKS**

Reduces risk of cross contamination of unused electrodes




# VersaTrode™

## Adult General Purpose Electrodes

Vermed VersaTrode™ electrodes by Nissha Medical Technologies are the ideal choice for adult cardiac monitoring and diagnostic applications. Designed to work in all procedures: stress test, long term monitoring, short term monitoring or resting EKG.

### A10005SG - 1.5" / 38.1mm - TEARDROP FOAM ELECTRODE

	Part #	Product #	Per Strip	Per Pouch	Per Case	Gel	Duration of Use
	6500021V	A10005SG-60	5	60	600		72 Hours
	6500484V	A10005SG-30	1	30	600		72 Hours
	6502954V	A10005SG-5NB	5	5	600		72 Hours
	6502947V	A10005SG-3NB	3	3	600		72 Hours
	6500015V	A10005-3-7%	3	3	600		72 Hours
	6500017V	A10005-7%NB	5	5	600		72 Hours
	6500019V	A10005-60-7%	5	60	600		72 Hours


### A10033 - 1.875" / 47.625mm - TEARDROP CLOTH ELECTRODE

	Part #	Product #	Per Strip	Per Pouch	Per Case	Gel	Duration of Use
	6500057V	A10033-1-60	1	60	600		72 Hours


### A10083 - 1.25" x 1.25" / 31.75MM x 31.75mm - SQUARE CLOTH ELECTRODE RADIOLUCENT & REPOSITIONABLE

	Part #	Product #	Per Strip	Per Pouch	Per Case	Gel	Duration of Use	MR Conditional
	6500075V	A10083-5	5	5	600		72 Hours	

### A10096 - 1.5" x 2.25" / 38.1mm x 57.15mm - RECTANGULAR FOAM ELECTRODE

	Part #	Product #	Per Strip	Per Pouch	Per Case	Gel	Duration of Use
	6502754V	A10096SG-30	1	30	750		72 Hours
	6502754V	A10096SG-3	1	30	750		72 Hours

### A10026T-60 - 1.875" / 47.625mm - TEARDROP TAPE ELECTRODE

	Part #	Product #	Per Strip	Per Pouch	Per Case	Gel	Duration of Use
	6500048V	A10026T-60	5	60	600		72 Hours