



ZBrand™ All Purpose Electrodes

FEATURES



VARIETY OF SHAPES AND SIZES
Available in multiple shapes and sizes for better nesting of electrodes in close spaces.



WET GEL OR SOLID GEL
Both gel options available to meet your clinical preference.



LIFT TABS
Lift tab allows for easy removal from liner.

ZBrand™ electrodes were designed for a general patient monitoring solution with quality performance.


PACKAGING

RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



ZBrand™

All Purpose Electrodes


The Vermed ZBrand™ Electrode line was developed to meet the requirements of the entire hospital by offering an economic, general patient monitoring solution with quality performance.

40 x 36 mm - Teardrop Cloth Electrode




Product #	Description	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502683V	Z05MC-3	3	3	1200	●		72 Hours
6502684V	Z05MC-5	5	5	1000	●		72 Hours
6502685V	Z05MC-50	5	50	1000	●		72 Hours

50 x 48 mm - Round Cloth Electrode



Product #	Description	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502686V	Z05MCRT-50	5	50	1000	●	⚠	72 Hours

40 x 36 mm - Teardrop Foam Electrode



Product #	Description	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502692V	Z05SG-3	3	3	1200	●		72 Hours
6502694V	Z05SG-5	5	5	1000	●		72 Hours
6502693V	Z05SG-30	3	30	1200	●		72 Hours
6502695V	Z05SG-50	5	50	1000	●		72 Hours
6502696V	Z05SGRT-30	3	30	1200	●	⚠	72 Hours
6502697V	Z05SGRT-50	5	50	1000	●	⚠	72 Hours
6502681V	Z05-5	5	5	1000	●		72 Hours
6502682V	Z05-50	5	50	1000	●		72 Hours

40 x 35 mm - Rectangle Foam Electrode



Product #	Description	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502721V	Z35-3	3	3	1200	●		72 Hours
6502722V	Z35-5	5	5	1000	●		72 Hours
6502723V	Z35-50	5	50	1000	●		72 Hours
6502724V	Z35RT-3	3	3	1200	●	⚠	72 Hours
6502725V	Z35RT-5	5	5	1000	●	⚠	72 Hours
6502726V	Z35RT-50	5	50	1000	●	⚠	72 Hours

40 x 36 mm - Teardrop Tape Electrode



Product #	Description	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502687V	Z05MT-3	3	3	1200	●		72 Hours
6502688V	Z05MT-5	5	5	1000	●		72 Hours
6502689V	Z05MT-50	5	50	1000	●		72 Hours
6502690V	Z05MTRT-5	5	5	1000	●	⚠	72 Hours
6502691V	Z05MTRT-50	5	50	1000	●	⚠	72 Hours