

Anderson Radiotransparent

Reference: D310

1 Pair / Pouch
10 Pairs / Case



PRODUCT INFORMATION

| | |
|---------------------------------|---|
| Pad Size | 4.69 x 6.13 in / 119.126 x 155.702 mm |
| Product Surface Area | 28.7497 in ² / 185.4815645 cm ² |
| Adhesive Area | 28.7497 in ² / 185.4816 cm ² |
| Conductive Area (2 Pads) | 32.1781 in ² / 207.6 cm ² |
| Gel Area | 16.6005 in ² / 107.1 cm ² |
| Pad Thickness | 0.094 in / 2.3876 mm |
| Lead Wire Length | 60 in / 1.524 m |

MATERIALS

| | |
|---------------------------------------|-------------------------------|
| Substrate Material | Polyethylene Foam |
| Adhesive | Acrylic, Medical Grade |
| Gel Characteristics | Conductive, Adhesive Hydrogel |
| Conductive Element | Carbon |
| Release Liner | Bleached Kraft Paper |
| Impedance Class | Low |
| Adaptor / Connector | Comparable to Philips OEM |
| Integrated Lead Wire Jacketing | PVC |
| Integrated Lead Wire Cord | Carbon |

ELECTRICAL PERFORMANCE (ANSI/AAMI EC 12)

Meets or Exceeds IEC 60601-2-4:2010

REGULATORY STATUS

| | |
|------------------|-----|
| CE Marked | Yes |
| 510K | Yes |



FEATURES

| | |
|---------------------------|------------------|
| MR Compatible | No |
| X-ray Translucence | Radiotransparent |

PACKAGING

| | |
|---------------------------------------|---|
| Product Weight | 2.4oz. / 68.0389g |
| Product Packaging Size (L x W) | 6.19 x 9.25 in 157.226 x 234.95 cm |
| Case Dimensions- Box (Size) | 10.75 x 6.75 x 6.5 in 273.05 x 171.45 x 165.1 cm |
| Leads Packaged | Inside Package |

BIOCOMPATIBILITY

| | |
|-------------------|--------|
| ISO 10993 | Passed |
| Latex Free | Yes |

ENVIRONMENTAL

| | |
|---|------|
| Halogenated Hydrocarbon Content (e.g. PVC) | None |
| Phthalate Derivatives Content (e.g. DEHP) | None |
| RoHS Compliant | Yes |
| REACH Compliant | Yes |

SHELF LIFE

| | |
|----------------------------|--|
| Product Shelf Life | 36 Months |
| Temperature Storage | Min -4; Max 104 °F / Min -20; Max +40 °C |

Vermed Part Number: 6500760V